

Physiotechnica



“AQUASPA” CONTACTLESS HYDROMASSAGE COUCH



MADE IN RUSSIA



“AQUASPA” contactless hydromassage couch

Overall dimensions (couch)	no more 2450x1100x590 mm
Weight of the couch (without water)	no more 180 kg
Power consumption*	3,5 kW

(*) - model with electronic individual programming of massage



“AQUASPA” CONTACTLESS HYDROMASSAGE COUCH (BATH) IS THE LATEST RUSSIAN MEDICAL HIGH-TECH EQUIPMENT, ALLOWING TO USE ALL THE THERAPEUTIC EFFECTS OF CLASSICAL HYDROMASSAGE WITHOUT DIRECT CONTACT OF THE PATIENT WITH WATER.

Advantages of the new model of 2019:

- ▶ Optimal length and height of the couch allow comfortable accommodation on it for the patient of almost any height and physique
- ▶ Two types of bracket with an LCD monitor are designed for comfortable management of the procedure by the staff or the patient himself
- ▶ The couch comes with three types of membrane: a translucent membrane with a chromotherapy system, a rubber membrane of blue or black color. Membranes may withstand distributed load up to 150 kg.
- ▶ Improved system of hydronozzles has a more powerful massage effect both in area and in strength
- ▶ The updated software with a user-friendly interface is optimal for managing the procedure
- ▶ Efficient cooling system allows you to conduct procedures without fear of overheating.
- ▶ The couch can be made in any color version according to the RAL palette, as well as an airbrush drawing can be applied on an individual order
- ▶ Various equipment sets of “Aquaspa” couch are possible



“AQUASPA” contactless hydromassage couch is equipped with a water heating system in the range of 30-35°C, which contributes to complete relaxation of the patient, and in case of use of the couch for carrying out different types of wraps, soft heat effect contributes to prolongation and enhancement of the therapeutic effect of the procedure.

In terms of effective area of impact on the human body among such products, the couch “Aquaspa” has no analogues in the Russian market!



“Aquaspa” couch massage programs:

STANDARD

Sequential activation of each of the 6 massage zones with an interval of 23 sec. Hydromassage begins with the upper shoulder girdle zone, successive movement of massage jets along the whole body of the person through all possible massage points of each of the massage zones to the calf muscle zone.

WAVE

Sequential activation of each of the 6 massage zones with an interval of 16 sec. Hydromassage begins with the upper shoulder girdle zone, successive movement of massage jets along the whole body of the person through all possible massage points to the calf muscle zone.

LYMPHATIC DRAINAGE

involves the passage of 6 massage zones along the entire human body from the muscle area of the lower limbs to the area of the upper shoulder girdle through all exposable massage points (an reverse algorithm to the “Wave” mode).

INTENSIVE

pairwise activation of the zones: the muscles of the brachial and thoracic area, lumbar and sacral area and the muscle area of the lower extremities with an interval of 23 seconds.

PULSE

Sequential activation of each of the 6 massage zones with an interval of 15 seconds, when moving from one zone to the next, the previous zone is activated 3 times for 1 sec in 1 sec. Hydromassage begins with the upper shoulder girdle area, successive movement of massage jets along the whole body of the person through all possible massage points.

FITNESS

a combination of the above modes within a set period of time.

INDIVIDUAL

involves selection of simultaneous operation of any one or two of the 6 massage zones within the specified time by the operator.



Application in medical practice:

the device is successfully used in health resort and rehabilitation institutions, spa salons and recreation complexes, centers of aesthetic medicine and cosmetology, in rooms of psycho-emotional relief and relaxation, rehabilitation centers for the disabled, in the sports medicine for relaxation after intense exercise.

Use effects:

- ▶▶ spasmed muscles are relaxed, thanks to it discomfort and pain disappear;
- ▶▶ blood supply of massaged tissues (skin, subcutaneous tissue, muscle and joint system) is improved, accordingly, their nutrition and oxygen supply are ameliorated;
- ▶▶ the metabolic processes in the body are activated, which has a beneficial effect on the overall condition of the body, the consumption of calories increases significantly. This is especially useful in weight loss programs;
- ▶▶ venous and lymphatic outflow is restored, which helps to eliminate edema and remove toxins from the body. This property is indispensable in chronic inflammatory diseases of the musculoskeletal system, as well as in programs of purification and detoxification;
- ▶▶ the function of internal organs is normalized through implementation of reflex mechanisms when paravertebral zones and biologically active points are affected - ideal for complex treatment of chronic diseases of cardiovascular, respiratory, digestive and other systems;
- ▶▶ the relief of fatigue and tension, restoration of psycho-emotional balance contribute to improvement of mood and general well-being, increase of vitality, rapid recovery of the working capability.

Aquaspa

“AQUASPA” contactless hydromassage couch (with electronic individual programming of massage)

“AQUASPA” contactless hydromassage couch (with the function of simultaneous hydromassage at 12 points)

“AQUASPA” heated water couch for therapeutic wraps

Technical specifications:

Overall dimensions (total) LxWxH, mm	2450x1100x1160	2450x1100x590	2450x1100x590
Overall dimensions (couch) LxWxH, mm	2450x1100x590	2450x1100x590	2450x1100x590
Weight of the couch (without water), kg	not more than 180	not more than 180	not more than 180
Total amount of water to be filled, l	220±10%	220±10%	220±10%
Working temperature of water poured into the tank of the couch, °C	30-35	30-35	30-35
Water heating time from room temperature (18°C) to 32°C, min	120	120	120
Power consumption, kW	3,5	3,0	2
Number of hydromassage programs	7	1	-
Number of hydronozzles	24	12	-
Number of hydromassage zones	6	6	-
Maximum distributed load on the bed up to, kg	150	150	150
Rated mains voltage at 50 Hz, V	220±10%	220±10%	220±10%
Programmable zone massage	+	-	-
Programmable bed heating mode	+	-	-
LCD touch screen	+	+	-
Blue, black 3 mm-thick bed/translucent bed with chromotherapy system	+	+	+
Possible additional options			
Single-step anti-slip stair-step	+	+	+
Indication			
Thalassotherapy wraps	+	+	+
Contactless massage	+	+	-

Physiotecnica

Russia, Saint Petersburg
 197198, Syezhinskaya Str. 23, lit. A, room 2-H
 Sales department: tel.: +7 (812) 321-67-80; 8-800-511-67-80
 e-mail: mail@pt-med.ru
www.pt-med.ru/en/

The booklet information is for informational purposes only